



- ★ Continue to care for the needs of your baby, including feeding and diaper changes.
- ★ Always provide safe sleep positions for your baby. Babies sleep best ALONE, flat on their BACKS, and in a CRIB without anything in the sleep area.
- ★ Protect your baby's eyes from the phototherapy using the provided eye protection.
- ★ Regularly check on your baby and monitor your baby's temperature.
- ★ Document your baby's temperature every 2-4 hours on the provided RECORD SHEET.
- ★ Do not put any creams or lotions on your baby's skin before or during phototherapy treatment.
- ★ Keep your baby flat on the nest mattress which provides the best possible treatment.
- ★ The infant's CRIB or PACK 'N PLAY is the best place for the bili-hut™ which prevents animals, siblings, or others from damaging the device which could cause harm.
- ★ The Environment High Temperature Warning automatically turns on and is **YELLOW** when the environment temperature surrounding the bili-hut™ is too high. The **BLUE** treatment lights turn off and the **YELLOW** indicator light turns on.  
(See additional instructions provided if Environment High Temperature Warning turns on)

### Important information regarding the bili-hut™

- ★ Never lay anything on top of the assembled bili-hut™ or its components.
- ★ Do not block the front and/or rear openings during treatment.
- ★ Do not use or store bili-hut™ in direct sunlight. This will cause it to overheat.
- ★ Do not transport the bili-hut™ with a baby inside.
- ★ Do not place anything inside the bili-hut™, including diapers, wipes, blankets, etc.
- ★ Do not place the bili-hut™ near a space heater, heating vents, or other device that gives off heat or has a hot surface.
- ★ Do not place the bili-hut™ near cooling sources such as fans, air conditioning or other devices that are cold or cool the air.
- ★ Do not place the bili-hut™ on an uneven surface or smaller than 2x4 feet surface without barriers.
- ★ Do Not place the bili-hut™ near a source of liquid that could spill into or onto the device.

Office Number 205-654-5339 (this number does not receive texts)