



Follow these Steps if the Environment High Temperature Warning Turns ON (**YELLOW**)



REMOVE YOUR BABY:

1. DO NOT turn off the bili-hut™.
2. Slide the nest out of the assembled hut.
3. Carefully lift your baby out of the nest and place them in a safe location.
4. Remove baby's eye protection.
5. Take your baby's temperature and document it on the Record Sheet.



1. Move the bili-hut™ away from any heat source such as heating vents.
2. Ensure the bili-hut™ is not in direct sunlight.
3. Lower the temperature in the room.
4. Move the bili-hut™ to a cooler, safe location.



Wait for the **BLUE** treatment lights to turn back on. This may take a few minutes. The bili-hut™ must be left ON during this time to see when the Environment High Temperature Warning automatically turns off.



The **YELLOW** indicator light turns OFF and the **BLUE** treatment light turns back ON.

RESUME TREATMENT

Place eye protection over your baby's eyes. Place your baby in the nest and gently slide the nest into the assembled hut until the nest touches the bumper at the foot-end of the base mattress.



If the **BLUE** treatment light remains off for >1 hour and you completed all possible troubleshooting items listed above call Enlighten Infant Home Care at 205-654-5339